

# ***Fish Tank Disinfection & Sterilization Algae Removal***

## ***Model: Doctor Aquatic Studio***



### **Features:**

#### **1. Inhibit green algae and promote plants growth**

- Inhibit phytoplankton initial growth and reproduction.
- Activate microelements to promote plants growth.
- Promotion of plants growth activates photosynthesis to stabilize ecosystem.

#### **2. Kill the crustacean pathogenic bacteria**

- Automatically kill the pathogenic bacteria in the water in each hour through direct and indirect disinfection, which greatly improves the survival rate of juvenile prawns.

#### **3. Depend on strong sterilization capabilities effectively kill pathogens**

- Aeromonas hydrophila, herpes, fungal, coli forms and other 99.9% kill efficiency results test.

#### **4. Safe because of coexistence with filter bacteria**



- As there is no residual toxicity, it is limited to aerobic bacteria.
- Compared to UV sterilizer, Chihiros DOCTOR does not work steadily so regenerating and recovery of filter bacteria is rapid.

#### **5. ALC program is intelligent management algorithm**


- Aquarium environment changes according to initial, middle, last stage, sterilizing intensity and operation cycle is intelligently designed to change by time series.
- As ALC program can simply choose a mode according to an aquarium capacity, it is very easy to use.

#### **6. The new Chihiros Doctor has integrated 3 application covering usage of aquatic plants, fish and shrimp keeping into only a single Doctor device.**

## Product Details:



## Installation & Connection:



## Operation Instructions:

①



- ◆ Press the power switch for long to start up(it is default to enter the function used for aquatic plants after starting up)
- ◆ Press the power switch for long to conduct function selection(used for aquatic plants,fish and crystal shrimps)

②



- ◆ Touch the power switch gently to select relevant capacity based on the capacity of the aquarium(if the capacity of the aquarium is just at the critical value,the small capacity is suggested to be chosen)

*Made in China*